

Laboratory Freezers B Medical Systems | F Range

Laboratory Freezers are devices intended for the storage of laboratory products, samples, test materials, chemicals or reagents at temperatures below -20°C.

Compliant to DIN 13221 | Medical Device Class I according to MDD 93/42/EEC



medical systems

SAVING LIVES THROUGH RELIABLE AND INNOVATIVE TECHNOLOGY

Medical Refrigeration



Safety Standards B Medical Systems



The Safety Standards developed by B Medical Systems define certain significant technical features of a product. These ensure the safe storage of the preparations as well as setting the highest standards of safety for the user.

PRECISION LINE	PREMIUM LINE	
		B Medical Systems Electronics
		B Medical Systems digital display
		B Medical Systems 4,3" digital color display with capacitive buttons
	\bigcirc	B Medical Systems 7" full touchscreen display
		Safety door lock (with 2 keys) and key-operated power switch ON/OFF (with 2 keys)
		Power indicator light and digital temperature indicator (display: 0.1 digits)
		Controlled fan cooling system for constant temperature and even temperature distribution across the entire refrigerating chamber
		Automatic switch-off of the evaporator fans when the door opens
		Self-contained alarm system with integrated battery takes over the alarm function and temperature value measurements in case of system failure for at least 48 hours
		Acoustic and visual alarm signal in case of temperature alarm and system failure
		The alarm history function on the electronic stores all the relevant values during a temperature alarm, such as: min., max., average temperature and alarm duration
		Remote transmission alarm signal (via potential-free contact) in case of temperature alarm (change-over contact)
		Door opening alarm (visual / acoustic)
		Designed and tested for climatic class SN (ambient temperature range +10°C to +32°C)
		Interior made from stainless steel
		Additional remote transmission alarm signal (via potential-free contact) in case of system failure (change-over contact)
		Smooth castors for optimum flexibility of movement
	\bigcirc	External water cooling
	\bigcirc	Ambient temperature sensor
0		B Connected - Universal software for the monitoring of refrigeration devices, including the acquisition, recording and visualization of temperature data
0		DCU (Data Communication Unit) - Hardware module monitoring all operating conditions and passing them through to a central database – via local network (in combination with °B Connected)
		Ethernet interface for the visualization of all operating and control functions (hardware and software settings) via °B Connected monitoring software on a peripheral device (computer)



Laboratory Freezers B Medical Systems | F Range

Laboratory Freezers are devices intended for the storage of laboratory products, samples, test materials, chemicals or reagents at temperatures below -20°C developed as a result of over 35 years of expertise in refrigeration technology.

7 models • Volume 121 > 895 L • Set temperature -41°C / -32°C • Climate class SN • Compliant to DIN 13221 | MDD 93/42/EEC, Class I

In conformity with national and international guidelines, regulations for Medical Devices offering reliability, efficiency and safety at an optimal price.



MODELS F130-290-380 | PRECISION LINE

• Integrated multi-function monitoring electronics

Integrated B Medical Systems multi-function electronics; with easy access in the door level; offers alarms, central alarm system, set point security, battery backup, full display and many other control features; compatible with °B Connected monitored solution for data monitoring and recording.

• Lifetime warranty corrosion-free cabinet

 One piece material cabinet with a lifetime warranty for zero corrosion, highest quality polypropylene for better insulation and no sharp edges to provide easier maintenance and cleaning

 Smart ergonomy - Easier and safer handling because of heavier components placed at the bottom

• High storage capacity

High storage per square meter in terms of net volume, with a very versatile inner volume utility due to multiple user friendly rails allowing modularity.



• Exclusive integrated electronics handle bar (Patent pending)

All functionalities are easily accessible:

- Digital control with full functionalities -Standard multifunction digital sensitive 4.3" panel integrated at optimal level in the door handle with full commands and connections, or optional full 7" touchscreen with pre-installed connection allowing exclusive °B Connected monitoring functionalities
- Open/close Heavy duty door lock mechanism

• High quality materials

New high quality coating, certified medical devices quality and antibacterial, high quality steel for better longevity and easy hygiene control.

• Versatile modularity

Excellent storage capacity and modularity - Large choice of shelves and drawers offering modularity for every need.

• Special features for easy deicing and defrosting

- New evaporator offering better performance than market average, resulting in lesser ice formation with automatic defrost
- Insulated inner doors for significantly lower loss of cooling when open



F RANGE Technical Data General features	F130	F290	F380	F400	F500	F700	F900
Gross / Net volume (I)	121/106	297 / 269	362 / 329	451/384	598 / 510	747 / 636	895 / 763
Gross / Net volume (I) Set temperature (preset)	121 / 106	297 / 269 -32°C	362 / 329	451 / 384	598 / 510	747 / 636 -41°C	895 / 763
	121 / 106		362 / 329	451 / 384			895 / 763
Set temperature (preset) Set temperature (setting range)	121/106	-32°C	362 / 329	451/384	-41	-41°C	895 / 763
Set temperature (preset) Set temperature (setting range) can be adjusted in steps of 0.1°C	121 / 106 42 min.	-32°C -32°C to -20°C	362 / 329 36 min.	451/384	-41	-41°C °C to -20°C	895 / 763

Climate class (ambient temperature range)		SN (+10°C to +32°C)
Defrosting technique	Manual	Automatic (hot gas evaporator defrosting)
Refrigerant type		R290

External dimensions H x W x D (mm)	830 x 595 x 695	1670 x 595 x 695	1980 x 595 x 695	1988 x 699 x 1039	1988 x 845 x 1039	1988 x 992 x 1039	1988 x 1139 x 1039
Inner dimensions H x W x D (mm)	630 x 475 x 470	1470 x 475 x 470	1790 x 475 x 470	1167 x 447 x 713	1167 x 590 x 713	1167 x 740 x 713	1167 x 887 x 713
Net weight with standard equipment (kg)	58	96	107				

	Supply voltage (V)	220-240	115-127	220-240	115-127	220-240	115-127	230	220	115-127	230	220	115-127	230	220	115-127	230	220	115-127
ner.	Frequency (Hz)	50 / 60	60	50 / 60	60	50 / 60	60	50	60	60	50	60	60	50	60	60	50	60	60
v Sol	Power (W)	300	360	300	360	300	360												
10 LU	Energy consumption (kWh/24h)	1.35 / 1.28	1.98	2.19 / 1.92	2.04	2.63 / 2.40	2.59												
L.	Heat emission (Kcal/h)	53	65	78	58	94	89												
	Compressor running time (%)	37	36	54	32	65	54												
	Noise level (dB(A)) (at 1m height & 1m distance)	45 / 47	47	42 / 43	45	46 / 47	47												

RANG RANG Technic Specific	cal Data			
		F130	F290-380	F400-500-700-900
Regulation se	NSOR 0°C to +180°C) ± 0.3°C		NTC	PT1000
Display sensor Precision (from -80 in reference body , corning 200-5CST	0°C to +180°C) ± 0.3°C / 100 ml DOW		NTC	PT1000
Accu data / fu panel when sy	Inction time of the control vstem failure	12 V - 2	.3 AH / 48h	12 V - 7 AH / 48h
Relative humi	idity at +32°C		S	75%
Door insulatio	on (polyurethane)	80 1	nm PU	65-82 mm PU
Casing insulat	tion (polyurethane)	60 1	mm PU	77-82 mm PU
	Inner body	Polypropylene	e (copolymer UV)	Stainless steel (1.4301)
	Outer body	Polypropylene	e (copolymer UV)	DX51D + Z100 coated
	Inner door	Polypropylene (copolymer UV)	Aluminium	Stainless steel (1.4301)
Material	Outer door	Polypropylene	e (copolymer UV)	DX51D + Z100 coated
	Separate interior doors		-	Anodized aluminium
	Drawer	Polycarbona	ite, transparent	· ·
	Shelf		-	Wire DIN177, PA11 coated
European Me	dical Device Directive		93 / 42 / EEC, Class	S (models: 220-240 V - 50/60 Hz)
EMC directive	•		2014 / 30 / EU ((models: 220-240 V - 50/60 Hz)
Low voltage a	lirective		2014 / 35 / EU ((models: 220-240 V - 50/60 Hz)
FDA regulatio	n		Class I medical devices acc. to 21	CFR Part 862.2050 (models: 115-127 V - 60 Hz)
cCSAus certifi	cate		Yes (models	5: 115-127 V - 60 Hz)
GMP - Clean r	room classification		Class ISO	6 / EC GMP B

Subject to change without prior notice.

SAVING LIVES THROUGH RELIABLE AND INNOVATIVE TECHNOLOGY



Equipment Standard & optional

F

RANGE

		F130	F290	F380	F400	F500	F700	F900
B Medical	Digital display	•	•	•	-	-	-	-
Systems	4.3" digital color display	-	-	-	•	•	•	•
Electronics	7" full touchscreen display	/ -	-	-	0	0	0	0
	Drawer with Front Cover	• 2	6	• 7	-	-	-	-
	Drawer without Front Cover	-	-	• 1	-	-	-	-
Interior equipment	Bottom Drawer	• 1	• 1	• 1	-	-	-	-
	Wire Shelf	-	-	-	• 5	• 5	• 5	• 5
	ST-Shelf	-	-	-	\bigcirc	\bigcirc	\bigcirc	\bigcirc
Separate inter in order to minimiz	rior doors the loss of cold air	-	-	-	•	•	•	•
°B Connected	- Monitoring Software	0	0	0	•	•	•	•
DCU - Data Co 'in combination wi	pmmunication Unit th °B Connected)	0	0	0	-	-	-	-
Ethernet inter	face	-	-	-	•	•	•	•
Ambient temp	perature sensor	-	-	-	\bigcirc	\bigcirc	\bigcirc	\bigcirc
Potential-free in case of system fo	alarm contact ^{ailure}	•	•	•	•	•	•	•
Integrated pol	rt for external sensor ^{mer)}	•	•	•	•	•	•	•
External wate	r cooling	-	-	-	0	0	0	0
Rollers Smoo	oth castors with stabilizers	• 2 -	• 2 -	• 2 -	- • 4	- ● 4	- • 4	- • 4
Door hinge Rig	ght Left							





Accessories Standard & optional

	F130	F290	F380	F400	F500	F700	F900
Dividers Set (H 35 mm or H 100 mm) to facilitate the handling and storage management	<u> </u>	()≤6	()≤7	-	-	-	-
S-Rack (H 100 x W 113 x D 638 mm)	-	-	-	0	0	0	0
Temperature chart recorder pre mounted (in the form of a circular chart recorder)	\bigcirc	\bigcirc	\bigcirc	\bigcirc	\bigcirc	0	0



Full functionalities at a glance

B Medical Systems | Electronics

PRECISION LINE | Multifunctional electronics with digital display and easy access, in the door handle; compatible with °B Connected monitored solution.



• Key power switch O = OFF / I = ON

• To switch the unit on, turn the key switch to position "I". The green power LED lights up

-30.5°C -32.0°C


-33.0°C

O Digital display

• After the self-check, the temperature inside the device is displayed

During normal operation the display shows:

temperature displayed in °C (Celsius) or



O Function buttons

• UP / DOWN buttons: navigation within the parameters and modification of their values

• ENTER button: use the ENTER-key to confirm the actual selection. This applies to the parameter value. The change of a parameter only becomes effective when confirmed

• MUTE / BACK button: in case of an alarm, the acoustic alarm can be muted for a specified time. In Parameter-mode the same key acts a "BACK-button" to go one step back

O USB-port

A USB port is available for data export only. Measured values of the display sensor are stored for at least 30 days (ring buffer) together with a time stamp.



/**+**_/

~

-

II X



• Acoustic and visual alarm signal

When an alarm occurs, the red alarm LED lights up and an acoustic signal sounds. The corresponding alarm message appears on the display alternately with the measured inside temperature. If several alarms occur simultaneously, the messages will be shown alternately.

The new B Medical Systems Electronics

offers a wide range of adjustment and diagnostic facilities as well as additional protection / warning operations (via external alarm operations, histories and individual display signals).

PREMIUM LINE | Multi-function electronics with 4,3" digital color display and capacitive buttons, integrated at optimal level in the door handle with full commands and connections.



-41.0°C Regulation

• Power switch OFF / ON

• To switch on press the ON/OFF button for 1 second. The green LED will light up

 Switch off protected by password. The door will automatically unlock before powering down

O Digital color 4,3" display

During normal operation the display shows:

- At the top: current voltage, current time, defrosting and gasket heater
- Left-hand side: internal temperature displayed in °C (Celsius) or °F (Fahrenheit) in increments of 0.1°

 Right-hand side: temperature of the regulation sensor. Displaying this temperature can be disabled through the settings



OUSB-port & SD-card-slot

Discover the ultimate control experience with the **new B Medical Systems 7**" full touchscreen display

More comfortable, more functionalities, no compromise!

• Acoustic and visual alarm signal

When an alarm occurs, the red alarm LED lights up and an acoustic signal sounds. The corresponding alarm message appears on the display together with the measured inside temperature. If several alarms occur simultaneously, the messages will be shown alternately.

• Capacitive buttons

- ARROW buttons: select screen items with the arrow buttons. In parameter editing mode, the arrow buttons are used to change values
- ENTER button: activates the currently selected item
- BACK button: goes one level upwards in the menu tree structure. It also serves as "escape" button
- DOOR UNLOCK button: unlocks the door of the device
- MENU button: opens the menu screen
- ALARM MUTE button: deactivates the acoustic alarm for a predefined period of time



•

G

:=

II ×



-41.0°C 📑

Optimum control and protection

B Medical Systems | °B Connected & DCU

°B CONNECTED | MONITORING SOFTWARE

Universal software for the monitoring of refrigeration devices, including the acquisition, recording and visualization of temperature data.

- Unique monitoring software for the full range of Blood Management Solutions and Medical Refrigeration products
- Web-based interfaces for computers and mobile devices
- Modern design for simple and intuitive use
- Graphic display of temperature curves
- Integrated event and activity history of appliances' components
- Data recording on centralized database for long-term archiving
- Easy setting of specific alarm forwarding, via email or SMS
- Generation of reports compiling crucial data and activities
- Temperature and detailed device data export for third-party software
- Important price advantage compared to a traditional circular chart recorder and its spare parts
- REST API to access raw data directly from the database in read only mode

KEY BENEFITS:

- Unique interface for the monitoring of the full range of refrigeration devices
- Centralized database providing data access to entire customer network

	URATION	
B con	nected	+O USER MANAGEMENT
	ER MMI (HIGDNO)	 Password login Access right
ne :1 Device configura	ation	+
PREMIUM LINE		DEVICE SETTINGS
		 Device information Device & user configuration
A AMERAN A AMERAN A AMERAN A AMERAN A AMERAN A AMERAN		((ر٢٠))
CST ed. (m.e.)	88 mm 2	ALARM PROFILE
BATOR		 Parameters of the alarm profile (day, hour, minute) Alarm management

2 DISPLAYING




DETAILED DISPLAY

 Optimised design for mobile and touch compatible devices
 Support of multiple languages

 Real time data: Current temperatures, set point, upper and lower limit, door state, compressor runtime
 Log data: Events, alarms within date/time, event type, description of the event. Event and alarm confirmation function
 Batch management: Read in of plasma bags with a barcode scanner
 Movement state and door control of Agitator devices
 Parameter settings of each device: Factory settings, real time clock, language, firmware version. R5485 address. IP address °B Connected & DCU offer a highly flexible system that is adaptable to specific customer requirements

- > Complete & legally safe documentation of temperature data
- > Comprehensive applications and diagnostic possibilities



DCU DATA COMMUNICATION UNIT (PRECISION LINE)

Hardware module monitoring all operating conditions and passing them through to a central database – via local network.



- Interface connection of B Medical Systems appliances to an existing network
- Graphical user interface displaying temperature and configurating devices and alarms easily
- Direct connection to Ethernet and serial bus RS485
- Digital in/out (programmable and customer-specific use)
- Recording and storing appliance relevant data
- Integrated USB port allowing to export stored data
- The DCU combined with the °B Connected software replaces the paper temperature recorder
- All data are recorded and save in the internal storage of the DCU and backed up in °B Connected data base if connected
- Several additional self-sufficient temperature sensors (up to 4 PT1000) may be connected
- Humidity sensor input (4-20 mA)

KEY BENEFITS:

- Central system collecting relevant temperature data of the appliances and their respective operating conditions
- A number of connection facilities allowing flexible upgrades for individual projects



*

Safe global blood management: from collection to transfusion, transportation, processing and storage



Reliable solutions for safe vaccination around the world



State-of-the-art technology for the exacting needs of the medical world

Our Global Expertise



After Sales support and service

We strive to provide you with the highest standards of service; not only through our selected distributors and partners for all your maintenance and service but also our second line trouble shooting and after sales service. This Factory based group of engineers is there to help our partners and yourself to get the best solution for your cold storage needs.



SAVING LIVES THROUGH RELIABLE AND INNOVATIVE TECHNOLOGY

B Medical Systems (formerly Dometic Medical Systems) has more than 35 years' experience in the medical refrigeration sector.

The company, formerly known as Electrolux Medical Systems, was founded in 1979 when the World Health Organization approached Electrolux in Vianden, Luxembourg, to create a solution for the safe storage and transport of vaccines around the world. In 2001, Electrolux Medical Systems became part of the Dometic Group, and was renamed Dometic Medical Systems. Having established a legitimate reputation in the medical equipment industry, the company has also become a global leader in vaccine cold chain.

B Medical Systems S.à r.l.

17, op der Hei L - 9809 Hosingen, Luxembourg

Tel.: (+352) 92 07 31-1 Fax: (+352) 92 07 31-300 info@bmedicalsystems.com

without prior notice.









Ed 10/2018