

**Don't just seal, SURESEAL**



**CENTRON**

*SureSeal™* SE160

TUBE SEALER

### **Versatility comes in standard**

*SureSeal™* SE160 is designed in consideration of multi tube sealing environment. It provides an ideal space-efficient solution.

#### ■ Features

- Premium sealing quality - new wide sealing technology
- Multi tube sealer set up - automatic sequential sealing
- Variable sealing distance : 70mm - 90mm (5mm interval)
- Simple operation, fast sealing, no warm-up required
- Automatic tubing detection
- LED indicator to show sealing progress
- Applicable to wide range of blood bag tubing
- Head protector / Splash guard for user safety

**Centron**

Smart Solutions for Smarter Bloodbanking

## New state-of-the-art sealing head

- New sealing technology ensures safe sealing of even thick and sensitive tubings with wide sealed surface.
- With automatic tubing detection, sealing will automatically start when tubing is placed in the sealing head.
- Semitransparent head protector for safe and easy operation (detachable, easy to clean)
- Applicable tubing diameter 2-6mm (1/4")




### Splashguard

- protects users from accidental blood splash
- Detachable, easy to clean

### Indication lamp

- shows sealing progress




### Wide sealed surface with notch at the center

- Reliable sealing quality prevents leakage
- Easy to segment



## Multi tube sealing set up example

- Linkage bars are provided for easy multi sealing head unit installation.



## Specification


Product	TUBE SEALER
Model	SE160
Main Power	100-240V~ 50/60Hz 200W
Dimension	69(W)mm/350(L)mm/245(H)mm

## Information



PENDING

## DISTRIBUTOR



**Centron** Centron Corporation

**Address:** 44 Beolmal-ro, Dongan-gu, Anyang, Gyeonggi-do,  
Postal Code 431-804 / Republic of Korea  
**Phone:** +82-31-422-7890 / **Fax:** +82-31-422-7866  
**E-mail:** info@centrons.com  
**Website:** www.centrons.com



※ The specification and design are subject to change without prior notice for product upgrades.